

**new from gSource**

**flush wire cutter with  
tc and silicone inserts**

### **Designed with safety and ease in mind.**

The gSource Flush Wire Cutter with tungsten carbide and silicone inserts can help prevent a cut piece of wire from being projected into the air or falling into the wound site.

- Silicone inserts are designed to hold the remnant wire for safe disposal after cutting.
- Improved design of tungsten carbide (TC) jaws cuts wire flush to the bone.
- Maximum leverage is achieved with the combined double action and leaf spring design. Provides a smooth and easy cutting action.
- Cuts k-wire up to OD .062" (1.6mm).
- End and side cutting jaws.
- Grooved handles provide a secure grip.
- Silicone inserts are suitable for use in manufacturing of medical devices. They are autoclavable and replaceable.
- Made from German stainless steel.
- Available from stock.
- Quality verified and guaranteed by gSource.

**gS 83.8450** 7", max cap .062" (1.6mm)

#### **Flush Wire Cutter with TC and Silicone Inserts**

Available from your gSource Distributor:

