


Introducing our new Zielke ring curettes for a comfortable grip.

The Zielke ring curette is ergonomically designed to provide a comfortable, secure grip during surgical procedures.

- White plastic handle is lightweight and 4 1/2" in length.
- Shaft is 9" in length and made from German stainless steel.
- Oval fenestrated cups are designed to provide better viewing of the surgical site. Available in 8mm and 11mm sizes.
- Shafts available in straight, angled and angled left or right patterns.
- Overall length 13 1/2"
- Quality verified and guaranteed by gSource.



gS 51.5640	straight	8mm	sharp/blunt	
gS 51.5645	straight	11mm	sharp/blunt	
gS 51.5650	angled 35°	8mm	sharp/sharp	
gS 51.5655	angled 35°	11mm	sharp/sharp	
gS 51.5660	angled 35°	8mm	sharp/sharp	right curved shaft
gS 51.5665	angled 35°	11mm	sharp/sharp	right curved shaft
gS 51.5670	angled 35°	8mm	sharp/sharp	left curved shaft
gS 51.5675	angled 35°	11mm	sharp/sharp	left curved shaft

Zielke Ring Curette

Available from your gSource Distributor: