

**Organize, store, and quickly access different sizes and styles of screwdriver bits.**

The gRack holds 3 1/2" SQC (small quick coupling) star and hex bits and a 4 1/2" SQC handle.

- Compact design - overall size is: 5 3/4" (h) x 7 1/4" (w) x 1 1/4" (d).
- Rack folds to close for convenient storage.
- When open, converts to a table top stand for use in O.R.
- Sizes are marked inside rack for easy identification.
- Made from anodized aluminum which is lightweight, yet built to withstand rigorous and repeated use.
- gSource SQC star bits, SQC hex bits and handle are sold separately.


Accommodates 1 each of the following:

SQC Handle	OAL (Overall Length)
gS 86.0040 black plastic and stainless steel	4 1/2"


SQC Star Bits	Size	OAL
gS 86.1604	T4	3 1/2"
gS 86.1608	T8	3 1/2"
gS 86.1715	T15	3 1/2"
gS 86.1725	T25	3 1/2"

SQC Hex Bits	Diameter	OAL
gS 86.1915	1.5mm	3 1/2"
gS 86.1925	2.5mm	3 1/2"
gS 86.1935	3.5mm	3 1/2"
gS 86.1940	4.0mm	3 1/2"

**gS 98.4050** holds 4 SQC star bits, 4 SQC hex bits, 1 SQC handle (all sold separately)



Closed position for storage. Empty



Open position as a table top stand. Loaded

## **gRack, Screwdriver Bits SQC**

For more information contact:

**gSource, LLC**  
 19 Bland Street  
 Emerson, NJ 07630  
 USA

P (800) 978-1119  
 (201) 599-2277  
 F (201) 599-3306  
 E email@gSource.com  
 www.gSource.com



ISO 13485: 2003 Certified  
 FS 589741  
 FDA Registered