

News from gSource - March 2019

[View this email in your browser](#)



News from gSource | March 2019

gSource is committed to putting the finest instruments into the hands of surgeons and their teams.

New Instruments

[Pin Shears Double Action 11 1/2"](#) [gS 83.9140](#)

- Shears 3.5mm [9/64"] and 4.0mm [5/32"] diameter titanium alloy pins with ease and provides a clean cut without burrs or sharp edges.
- Grooved handles are ergonomically designed for a comfortable and secure grip.
- Double action design allows for ease in cutting.



Feature Instrument

Clean Wave Kerrison Punch Runner-Up in Surgical Products Award

In the November/December issue of Surgical Products magazine, the gSource Clean Wave Kerrison Punch (gS 69.3081-gS 69.3085) was announced as a runner-up for their Excellence in Surgical Products (ESP) award in the Instrumentation category.

- The wavelike shaped design allows cleaning utensils to easily reach through the recesses and openings between the main body and slider for cleaning prior to sterilization.
- Wavelike shape also helps to reduce friction between the slider and main body, providing smooth function.
- 8" shaft with forward angled foot plate.
- Forward angled foot plate provides precise and controlled cutting action.
- Features a thin foot plate design.
- Available in 1mm-5mm bite sizes.
- Bone ejector is incorporated into the design of 2mm-5mm bite sizes.
- Satin finish helps to reduce glare.



Services

Custom Instrumentation

gSource will create entirely new instruments in a timely and economical manner. For an accurate quotation and delivery schedule, send us a sample, sketch or drawing. Custom marking, packaging and labeling can also be provided.



Instrument Repair Services

gSource offers sharpening, spring and screw replacement and refurbishing repair services to ensure quality performance during the entire life of your gSource instruments.

Our in-house repair technicians will expertly sharpen and repair your gSource instruments, and other instrument brands, according to gSource Verified Quality standards. Our superior customer service teamed with our highly skilled instrument repair technicians will ensure you receive an excellent repair experience.

See www.gSource.com for detailed information describing our repair services and to print/view our [Order Form](#) and [Price List](#). If you have any questions about our repair services, please contact the gSource Repair Center at 800-978-1119 or repair@gSource.com.



Trade Shows

Exhibiting at:

AAOS 2019 Annual Meeting
March 13 – 15, Las Vegas, Nevada

AAOS
AMERICAN ACADEMY OF
ORTHOPAEDIC SURGEONS

gSource Booth 1322

OMTEC 2019
EXHIBITOR

June 12-13, 2019
Rosemont, IL

gSource Booth 135

