

Bend rods with ease and comfort!

- Bends rods up to 4.5mm with three different cam settings.
- Separate reduction ring not required.
- Adjustable center cam is spring loaded.
- Biocompatible grooved silicone inserts on handles help to prevent slippage and provide a secure grip.
- Silicone inserts have excellent resistance to alkaline, saline and acidic solutions.
- Non-fading silicone color withstands repeated autoclaving.
- Made from German stainless steel.
- Quality verified and guaranteed by gSource.

gS 82.7695 11", for rods up to 4.5mm [.177"]
green silicone inserts on handle

Universal Rod Bender with silicone inserts

For more information contact:

gSource, LLC
19 Bland Street
Emerson, NJ 07630
USA

P (800) 978-1119
(201) 599-2277
F (201) 599-3306
E email@gSource.com
www.gSource.com



Handle showing silicone insert.
Detail not to scale.

