

shearing action cutters



shear without burrs or sharp edges!

gSource[®]
An ARCH Medical Solutions Company

2 - realistic price

3 functions in 1 versatile instrument!

- 1) Shears wires with ease and provides a clean cut without burrs or sharp edges.
- 2) Bends wires quickly and easily.
- 3) Cross-serrated TC inserts ensure a secure grip on wires for pulling.


Features:

- Tungsten carbide (TC) in jaws and cutter.
- Grooved handles are ergonomically designed for a comfortable and secure grip.
- Made from German stainless steel.
- Quality verified and guaranteed by gSource.

gS 81.3380 8"


Wire Bending Pliers with Cutter

tungsten carbide jaw
max cap 2.0mm [5/64"]



Shear 3 different pin diameters with 1 instrument!

- Shears 2.5mm [.098"] 3.0mm [.118"] and 3.5mm [9/64"] diameter pins with ease and provides a clean cut without burrs or sharp edges.
- Grooved handles are ergonomically designed for a comfortable and secure grip.
- Double action design allows for ease in cutting.
- Made from German stainless steel.
- Quality verified and guaranteed by gSource.



gS 83.9135 11 1/2"


Pin Shears DA

for 2.5mm, 3.0mm, 3.5mm diameter pins



Shear 3 different rod diameters with 1 instrument!

- Shears 4.5mm [.177"], 5.5mm [.217"] and 6.35mm [1/4"] diameter rods, leaving a smooth and clean surface rather than a sharp, jagged, burr-like surface common when using a standard pinching-type rod cutter.
- Available with and without detachable handles.
- Double action design allows for ease in cutting.
- Biocompatible silicone handles help to prevent slippage and provide a secure and comfortable grip.
- Silicone handles have excellent resistance to alkaline, saline and acidic solutions.
- Seven standard biocompatible silicone colors are available or a custom color of your choice.
- Non-fading vibrant silicone colors withstand repeated autoclaving.
- Made from German stainless steel.
- Quality verified and guaranteed by gSource.
- Note: Instrument is not designed to be used inside the patient's body.



gS 83.7280
non-detachable handles

gS 83.7280 22", non-detachable handles, black silicone
gS 83.7290 22", with detachable handles, black silicone

gSilicone Rod Shears DA
 for 4.5mm, 5.5mm, 6.35mm diameter rods
 black silicone handles

gSource

gSource is committed to putting the finest instruments into the hands of surgeons and their teams.

Not all instruments are created equal. Adequate for the job is, in reality, inadequate. A better instrument reduces distractions and facilitates operations. A better instrument helps a surgeon perform at his/her best. A better instrument leads to better results. gSource provides better instruments.

Orthopedic Focus

We are more than a source for quality instrumentation. gSource is an advocate for the orthopedic community. We can be relied upon to provide superior instrumentation for spinal and orthopedic procedures, and do so in a time critical fashion.

gSource Instrument Attributes

Whether crafted from German surgical stainless steel, or machined from U.S. surgical stainless steel, our instruments are recognized by their finely finished surface. They are designed to perform with precise surgical function and are also conceived to be affordable. It is this combination that distinguishes the gSource brand.

gSource Verified Quality

Every gSource instrument must pass detailed Quality Assurance tests before it can be sold. Instruments are tested for quality of material, workmanship, function, pattern consistency and critical dimensions.

Instrument reSource

The gSource catalog contains 3,500 instruments primarily used in orthopedic and spinal procedures. To view our catalog or for a quick answer on availability of other patterns, refer to our website at www.gSource.com. Please inquire about the availability of any instrument not shown in the catalog.

Custom Instrumentation

We will create entirely new instruments in a timely and economical manner. For an accurate quotation and delivery schedule, send us a drawing, sketch or sample.



gSource, LLC
19 Bland Street
Emerson, NJ 07630
USA

P (800) 978-1119
(201) 599-2277
F (201) 599-3306
E email@gSource.com
www.gSource.com

Service First

Your inquiries will be answered quickly and accurately by knowledgeable professionals.

Trusted Supplier to World Leaders

Many world leaders and innovators in orthopedics and spine have found in gSource a trusted and reliable partner. We provide superior instrumentation in a time-critical fashion and at realistic prices.

Guarantee

All standard instruments are guaranteed for life against manufacturing defects of material and workmanship. Any instrument proving to be defective will be replaced or the purchase price refunded.

1. Tungsten carbide inserts are guaranteed for three years. Replaceable parts, other than springs, are guaranteed for one year. Replaceable springs are guaranteed for life against manufacturing defects of material and workmanship.
2. This guarantee is void if instruments are altered or not maintained or repaired properly or if they are not used for their intended surgical purpose.
3. Any unused instrument may be returned for full credit within 90 days of invoice date.

gSource Advantages

- On-time delivery
- Realistic price
- Verified quality
- German craftsmanship
- Precise U.S. machining
- Orthopedic and spinal focus
- Product development support
- CAD support
- Custom labeling, packaging and marking
- Large selection and inventory
- Adaptability/flexibility
- Superior customer service
- Full satisfaction guarantee



ISO 13485: 2003 Certified
FS 589741
FDA Registered



gSource, LLC Founded 1999