

Comfort on hand!

- Ergonomically designed to provide a comfortable grip during surgical procedures.
- Biocompatible silicone helps to prevent slippage and provide a secure grip.
- Silicone has excellent resistance to alkaline, saline and acidic solutions.
- Non-fading vibrant RAL 5015 color withstands repeated autoclaving.

Specifications:



- 51mm: Overall Length
- 66mm: Outside Diameter
- 33mm: Bottom Inside Diameter
- 44mm: Inside Height
- 9.5mm: Bottom Opening Diameter
- Quality verified and guaranteed by gSource.



Outside View



Inside View



Price: \$75.00 each

Z 18040402 2 1/2", sky blue

Silicone Ball Handle

For more information contact:

gSource, LLC
19 Bland Street
Emerson, NJ 07630
USA

P (800) 978-1119
(201) 599-2277
F (201) 599-3306
E email@gSource.com
www.gSource.com



ISO 13485:2016 Certified
FS 589741
FDA Registered